

# MIXBAGS

## Sterile bags



Mixbag full filter



Mixbag side filter

Reference	Designation	Nominal volume	Useful volume	Filter porosity	Dimensions	Unit per bag	Bag per box
MIXBAG Full filter	MIXW3002	400ml	400ml	<150µm	190x300mm	500	3
MIXBAG Side filter	MIXW3003	400ml	400ml	<150µm	190x300mm	500	3

- Non woven filter
- Food contact: Compliant with FDA regulations and EU Directive 2002/72/EC
- Manufactured following the ISO9001-2000
- All MIXBAGS are sterilised by irradiations. A certificate of irradiation is available upon request. Each box is delivered with a red sticker on it



Fully compatible with  
ProBlend/ProDilute



#### Contact Synbiosis

##### Europe:

Beacon House Nuffield Road  
Cambridge CB4 1TF UK  
Tel: +44 (0)1223 727125  
Fax: +44 (0)1223 727101  
Email: [sales@synbiosis.com](mailto:sales@synbiosis.com)  
[www.synbiosis.com](http://www.synbiosis.com)

##### USA:

5103 Pegasus Court Suite L  
Frederick MD 21704 USA  
Tel: 800-686-4451 301-662-2863  
Fax: 301-631-3977  
Email: [ussales@synbiosis.com](mailto:ussales@synbiosis.com)